

Case Study #1 – AT Category: Hearing 1 (SETT)

Student – Pre-K student

Environments – Regular education classroom, gym, and art class

Tasks – Listening to instruction in variety of environments and varying noise levels

Tools – Does the student already have access to tools? There is no classroom amplification system, and the student has just been fitted with hearing aids

Low-Tech

Mid-Tech

High-Tech

Visual Aids – Notes / Captions – Lip Reading – iPod - FM Listening system – Digital Listening System – Roger System



STEP 1: Based on S-E-T data, enter descriptors or functions needed by the student across the shaded top row - 1 descriptor per column

STEP 2: Enter promising tools in the shaded left column - 1 tool per row

STEP 3: Note whether each tool matches a descriptor by placing an “X” in each of the applicable white boxes

Descriptors	Cordless	Separate Transmitter and Receiver	Rechargeable	Uses Ear buds, headphones, or neckloop	Able to be taken to multiple environments	Pairs with telecoil in hearing aid or implant	Non-interfering frequency	Improves noise to sound ratio
Tools								
Comfort Contego	X	X	X	X	X	X	X	X
Phonak MyLink w/Roger Inspiro	X	X	X	X	X	X	X	X
Phonak Roger X receiver with Roger Inspiro or Roger Pen	X	X	X		X		X	x
Redcat or other classroom audio system	X	X	X					
Comfort Duett or Pocket Talker	X		X	X		X	X	